

Developing an Oxygen Detection Device for a Microfluidic Hypoxia Chamber
Experimental Protocol and Procedures
November 16th, 2012

Client: *Professor Brenda Ogle, PhD*

Advisor: *Professor Randolph Ashton, PhD*

Team: *Matthew Zanutelli*
Chelsea Bledsoe
Karl Kabarowski
Evan Lange

EXPERIMENTAL PROCEDURE FOR CREATION OF STANDARDIZED CURVE AND OPTIMIZATION OF THIN-FILM SENSOR

1. Prepare PtOEPK matrix at different volume ratios of matrix: PtOEPK
 - a. Test to find optimal ratio
 - b. Dissolve PS pellets in toluene solution (7% w/w ideal, test dramatically above and below)
 - i. Possible tests: 17% w/w and 1% w/w
 - ii. GOAL: Find optimal w/w percent of encapsulation matrix.
 - c. PtOEPK dye added to matrix
 - i. Test conditions of 0.1mg: 1mL, 1mg: 1mL, and 10mg: 1mL
 - ii. GOAL: Find optimal ratio of luminescent material in encapsulation matrix
2. Transfer ~30 μ L of PtOEPK matrix into each well of the 96-well plate
 - a. Need to test to determine effective volume needed to completely cover bottom of well
 - b. Films will need to be as thin as possible to optimize response time
 - i. If time, attempt to test films of varying thickness
 - c. Spin-state wells to pull all material to bottom of wells
3. Allow the PtOEPK matrix to stand in wells for 24hrs at room temperature in dark for gelation and drying
 - a. NOTE: may require longer time to dry
 - i. Need to determine time needed for films to gel
4. Vary concentration of oxygen in 96-well plates (percentage of oxygen, v/v) by manipulating the flow rate of oxygen and nitrogen flowing into the microwell containing the oxygen sensing film
 - a. Vary concentrations from 1% to 25% if possible
5. The quenching response is determined by measuring the luminescence intensity of films under no oxygen and a fully oxygenated solution
6. Plot quenching response versus volume ratio
7. Determine Stern-Volmer plots by measuring luminescence intensity under various gaseous oxygen or dissolved oxygen concentrations.

**At this point in the design photostability and leaching effects will not be tested as they have thoroughly be proved by other labs and previous work

HYPOTHESIS:

From previous research and findings, we believe that the thin-film sensor will work most efficiently with 7% w/w PS in toluene solution and a 1mg: 1mL PtOEPK in encapsulation matrix. Additionally, as oxygen concentration is increased intensity of the luminescent material should decline as in the presence of oxygen oxidation occurs with Pt- causing luminescent properties to decline.

EXPERIMENTAL DESIGN

Test will need to be performed multiple times to demonstrate high reproducibility of the thin-film sensors

- Test 1A:
 - 1% w/w PS in toluene solution
 - 0.1 mg: 1mL ratio of PtOEPK in encapsulation matrix
- Test 1B:
 - 1% w/w PS in toluene solution
 - 1 mg: 1mL ration of PtOEPK in encapsulation matrix
- Test 1C:
 - 1% w/w PS in toluene solution
 - 10 mg: 1mL ration of PtOEPK in encapsulation matrix
- Test 2A:
 - 7% w/w PS in toluene solution
 - 0.1 mg: 1mL ratio of PtOEPK in encapsulation matrix
- Test 2B:
 - 7% w/w PS in toluene solution
 - 1 mg: 1mL ration of PtOEPK in encapsulation matrix
- Test 2C:
 - 7% w/w PS in toluene solution
 - 10 mg: 1mL ration of PtOEPK in encapsulation matrix
- Test 3A:
 - 17% w/w PS in toluene solution
 - 0.1 mg: 1mL ratio of PtOEPK in encapsulation matrix
- Test 3B:
 - 17% w/w PS in toluene solution
 - 1 mg: 1mL ration of PtOEPK in encapsulation matrix
- Test 3C:
 - 17% w/w PS in toluene solution
 - 10 mg: 1mL ration of PtOEPK in encapsulation matrix

**For all test conditions, quenching rate, max and min intensity, and Stern-Volmer rates will be determined

**Due to time constraints, team will begin with Test 2, as it is hypothesized to be the optimal conditions.

ANALYSIS

- Optimal composition of thin-film sensor should provide highest sensitivity
 - Sensitivity can be expressed by quenching response and quenching ratio (I_0/I_{100} , where I_0 is luminescence intensity at 100% N₂ and I_{100} is luminescent intensity at 100% O₂)
- Films will need to be as thin as possible to optimize response time to excitation
- Test will need to be performed multiple times to demonstrate high reproducibility of the thin-film sensors