

BIOMIMETIC BONE MARROW MICROENVIRONMENT FOR hMSC CULTURE



Taylor Marohl, Maddie Meier, Veronica Porubsky, and Michelle Tong
 Client: Professor Wan-Ju Li | Advisor: Dr. Tracy Puccinelli
 Department of Biomedical Engineering, University of Wisconsin-Madison



ABSTRACT

Degenerative musculoskeletal diseases, specifically osteoarthritis, affect approximately 1 out of 3 people aged 65 and over [1]. Noninvasive treatment approaches only alleviate symptoms and do not address the degeneration of joints [2,3]. As an alternative, our client is interested in employing a regenerative medicine approach with human mesenchymal stem cells (hMSCs). Thus, a new culture system is needed to maintain hMSCs in a quiescent state. Last semester, rheometer testing characterized bone marrow stiffness (E) within a range of 61.7- 312.8 kPa. However, we will proceed with literature values due to their consistency across studies. This semester, we fabricated RGD-conjugated PEG hydrogels with Young's moduli spanning the human bone marrow stiffness range. We also developed a bioreactor that maintains [CO₂] and [O₂] within ±1.5% of 1% and 5%, respectively. Future work will fabricate fibronectin-conjugated PEG and 2 additional bioreactor chambers, ultimately allowing for multivariate testing.

SIGNIFICANCE

Musculoskeletal Diseases:

- Autoimmune, inflammatory, and degenerative diseases
 - ex: osteoarthritis affects 33.6% of individuals age 65+ [1]
- Potential treatment: regenerative medicine with human mesenchymal stem cells (hMSCs)

Bone Marrow and hMSCs:

- Niche of quiescent hMSCs; visco-elastic tissue [4]
- Complex environment allows maintenance of hMSCs in undifferentiated state
- Stiffness: 0.1-0.2 kPa [5] Oxygen tension: 1-5% [14]
- In vitro* culture: loss of potency [5]





Figure 1. MSC source.

PROBLEM STATEMENT

Current culture methods are variable in their ability to maintain hMSCs in a multipotent and quiescent state, and do not attempt to **recreate** the **physiological conditions** that **prevent differentiation**. For this reason, there is a need for a culture system that sustains cell quiescence by mimicking the bone marrow niche through **substrate stiffness** and **oxygen tension**.

INNOVATION



- Cultures living bone marrow
- Retains hematopoietic stem cells and progenitor cells for one week
- Drawbacks:
 - Bone marrow generation in mouse
 - No consideration of oxygen tension or mechanical properties

Figure 2. Bone marrow-on-a-chip [6].

BIOREACTOR

Aim: Develop bioreactor that maintains user-set oxygen tension.

Design:

- Hypoxia chamber housed inside cell culture incubator, which provides temperature regulation




Figure 3. Bioreactor block diagram illustrating system function.

- Solenoid valves open for 2 seconds when CO₂ or O₂ concentrations drop below user-set values, close for 7 seconds.





Figure 4. Bioreactor prototype, including circuit with solenoid valves, Arduino, and gas sensors.


Validation:

- Monitor [CO₂] and [O₂] for 10 minutes once stable
- Variation within ±1.5% met for both gases, fine-tune flows to improve




FUTURE WORK

Scaffold: Fabricate fibronectin-conjugated PEG gel that allows a chemically-defined environment.



Bioreactor: Fabrication of 2 additional chambers to allow higher throughput evaluation of hMSC behavior at various oxygen tensions.

Cell Culture Evaluation: Culture hMSCs on scaffolds of varying stiffnesses, and on tissue culture plastic in varying oxygen tensions.




Multivariate Testing: Cell culture evaluation over a range of biomaterial stiffnesses and oxygen tensions could reveal the ideal combination of factors for maintaining hMSC quiescence *in vitro*.

SCAFFOLD

Aim: Fabricate RGD-conjugated PEG hydrogels with Young's moduli spanning relevant human bone marrow stiffness range.

Design:



Gel Attempted	Primary Troubleshooting Issues			
	Hydrogel Formation	Peptide Lyophilization	Peptide Dialysis	Adhesion
3D Porous PEG	X			
2D PEG-RGD		X	X	X
2D PEG-GRGDS			X	X
2D PEG-gelatin	X			X

Validation:

- Cell adhesion test with hMSCs
- Mass-spectrometry of peptide aminolysis products (future)
- Rheometry for Young's Moduli of 2.5 - 7.5 w/v% PEG gels (future)
- Experiment for effect of stiffness on quiescence (future)

ACKNOWLEDGEMENTS

Dr. Tracy Puccinelli, Professor Wan-Ju Li, Dr. John Puccinelli, Dr. Joseph Towles, Tony Berger, Ronald McCabe, and Kevin Wang

REFERENCES

[1] Osteoarthritis (OA) | Arthritis | CDC. *Cdc.gov*, 2016. [Online]. Available: <http://www.cdc.gov/arthritis/basics/osteoarthritis.htm>. [Accessed: 09- Dec- 2016].
 [2] A. Gelber, "Osteoarthritis," *Annals of Internal Medicine*, vol. 161, no. 1, pp. ITC1-1, 2014.
 [3] E. Barnes and N. Edwards, "Treatment of Osteoarthritis," *Southern Medical Journal*, vol. 98, no. 2, pp. 205-209, 2004.
 [4] L. E. Jansen, N. P. Birch, J. D. Schifman, A. J. Crosby, and S. R. Peyton, "Mechanics of intact bone marrow," *Journal of the Mechanical Behavior of Biomedical Materials*, vol. 50, pp. 229-307, 2015.
 [5] J. P. Winer, P. A. Janney, M. E. McCormick, and M. Funaki, "Bone marrow-derived human mesenchymal stem cells become quiescent on soft substrates but remain responsive to chemical or mechanical stimuli," *Tissue Engineering Part A*, vol. 15, no. 1, pp. 147-154, 2009.
 [6] Y. Torisawa et al. "Bone marrow-on-a-chip replicated hematopoietic niche physiology *in vitro*," *Nature Methods*, vol. 11, no. 6, pp. 663-669, 2014.
 [7] N. L. Petrakis, "Temperature of human bone marrow," *Journal of Applied Physiology*, vol. 4, no. 7, pp. 549-553, 1952.
 [8] T. A. Metzger, J. M. Shudick, R. Seekell, Y. Zhu, and G. L. Niebur, "Rheological behavior of fresh bone marrow and the effects of storage," *Journal of Mechanical Behavior and Biomedical Materials*, vol. 40, pp. 307-313, 2014.
 [9] Holzwarth C, Vaegler M, Giescke F, Pfister SM, Handgreiner R, Kerst G, Müller J. "Low physiologic oxygen tensions reduce proliferation and differentiation of human multipotent mesenchymal stromal cells." *BMC Cell Biology* 11:11, 2010.
 [10] C. Zhang, S. Hekmatfar, A. Ramanathan, and N. W. Karuri, "PEGylated human plasma fibronectin is proteolytically stable, supports cell adhesion, cell migration, focal adhesion assembly, and fibronectin fibrillogenesis," *Biotechnology Progress*, vol. 29, no. 2, pp. 493-504, 2013.
 [11] Zhou, Yaxian, Michael Zimber, Huihua Yuan, Gail K. Naughton, Ryan Fernan, and Wan-Ju Li. "Effects of Human Fibroblast-Derived Extracellular Matrix on Mesenchymal Stem Cells." *Stem Cell Reviews and Reports* 12.5 (2016): 560-72. Web.
 [12] Okolicsanyi, Rachel K., Emily T. Camilleri, Lotta E. Olkari, Chieh Yu, Simon M. Cool, Andre J. Van Wijnen, Lyn R. Griffiths, and Larisa M. Haupt. "Human Mesenchymal Stem Cells Retain Multilineage Differentiation Capacity Including Neural Marker Expression after Extended In Vitro Expansion." *PLoS ONE* 10.9 (2015): n. pag. Web.
 [13] ThermoFisher Scientific. (2016). *Heracell VIOS Incubator*. Retrieved from: <https://www.thermoFisher.com/order/catalog/product/51030284>
 [14] Mohtayidin, Ahmed et al. "Oxygen in Stem Cell Biology: A Critical Component of the Stem Cell Niche." *Cell Stem Cell* 7:2 pp. 150 - 161, 2010