AUTOMATED URETERO-INTESTINAL ANASTAMOSIS WITH ABSORBABLE STAPLES

John Cheadle, Rebecca Clayman, Terra Gahlman, Katie Pollock and Kim Safarik


Client: Dr. Tracy Downs, Department of Urology, Madison Advisor: Dr. Wan-Ju Li, Department of Orthopedics, UW Madison

March 4, 2011




BLADDER CANCER

- 5th most common cancer in US
- 70,530 new cases and 14,680 deaths in 2010
- Radical cystectomy
 - Required when muscle layer is invaded
- Urine storage
 - Urostomy bag or neobladder



CURRENT TREATMENT: BLADDER CANCER



http://www.ucurology.urologydomain.com/images/uploaded/whu/ileal1.gif

- Bladder removed
- Neobladder formed from small bowel
- Two ureters attached via suture

PROBLEM STATEMENT

- Ureters connected to new bladder tissue via suture
- Lengthy procedure time
- Metal staples = stones

• Goals


- Semester 1: Small stapler
- Semester 2: Continued stapler work with testing and staple research



http://www.swicofil.com/images/suture_needle.jpg


DESIGN REQUIREMENTS

- Open surgery
- Sterile
- 0.7cm diameter
- Operate with single motion
- Must create a water tight seal
- Mitigate tissue damage




http://bostinnovation.com/wp-content/uploads/2010/06/staple-icon.jpg


STAPLER PARTS : ANVIL AND NEEDLE



STAPLER PARTS: RING CLAMP



STAPLER PARTS: TEETH AND CARTRIDGE



STAPLER PARTS: FIRING MECHANISM

- 1. Two compressive springs
- 2. Washers
- **3**. Lever and pin
- 4. Flange translator
- **5**. Two concentric rods



STAPLE RESEARCH: COMPOSITION

- Ideal properties: rigid, but flexible; biodegradable
- Polylactic acid (PLA): imparts toughness and strength
- Polycaprolactone (PCL): imparts flexibility

Mixing ratios:

70/30 PLA:PCL with dicumyl peroxide (DCM)¹

- Exhibits plastic properties similar to ABS or HIPS¹
- Experiment with other ratios

STAPLE RESEARCH: PROCESSING

- Creation of PCL/ PLA polymer with DCP¹
 - Twin-screw extrusion





Fig 1. Top cross section

- Injection Molding
 - Forms polymer into staple shape for testing



STAPLE RESEARCH: SECURING TISSUE

Staple

 Bending testing
Ultrasound polymerization of ends




Fig 1. Staple being bent in ring clamp (cross section)

TESTING

Staples

- Testing multiple formulations
 - Qualitative do they bend?
 - Tensile testing
 - Compare with sutures
 - Polymerization testing

Stapler


- Prototype efficacy with metal staples
 - Synthetic material
 - Tissue explants
- Leak Testing

FUTURE WORK

- Testing with Absorbable Staples
- Publication in Journal
 - Journal of Endourology
- Application into radical cystectomy procedure



Ralph V. Clayman, M.D. Arthur D. Smith, M.D.



The Official Journal of the Endourological Society Society for Engineering and Urology Society of Urologie Robotic Surgery

Mary Ann Liebert, Inc. Epublishers

REFERENCES

• Semba, T, K Kitagawa, US Ishiaku, and H Hamada. "The Effect of Crosslinking on the Mechanical Properties of Polylactic Acid/polycaprolactone Blends." *Journal of Applied Polymer Science*, 101.3 (2006): 1816-1825.

ACKNOWLEDGEMENTS

- Dr. Tom Turng
- Professor Kim Manner
- Dr. Tim Oswald
- Dr. Tracy Downs
- o Dr. Wan-Ju Li

QUESTIONS?

