Non-Invasive Cervical Cancer Screening



Georgia Hancock Mira Baichoo Karina Buttram Josephine Hall Adrienne Simpson Cora Williams Team Leader Communicator BPAG BSAC BWIG BWIG

Georgia Hancock Human Papillomavirus (HPV) and Cervical Cancer

- HPV is the most common sexually transmitted infection
 - 200 different types, but only 40 can infect the genital areas
- Persisting infection of certain strains of HPV can lead to cervical cancer
 - HPV 16 is linked to approximately 50% of cervical cancers worldwide
 - HPV 18 is the second most important
- Cervical cancer is one of the most common cancers in women





Figure 1: Computer generated model of HPV

Georgia Hancock Current Cervical Cancer Screening Procedures



- Provider inserts speculum to visualize cervix
- Provider then uses a wooden or plastic scraper and/or cervical brush to collect cell sample
- A cytopathologist then examines the cell sample under a microscope

@ MAYO FOUNDATION FOR MEDICAL EDUCATION AND RESEARCH. ALL RIGHTS RESERVED.

Figure 2: Step-by-step illustration of a Pap smear test

[•] Pap smears

Cora Williams

Problem Statement

- Cervical cancer is one of the most treatable cancers when diagnosed early
- Current screening methods must be performed by medical professional
 - Not easily accessible in developing countries
 - Uncomfortable experience
- Discrete, self-collected urine sample test
 - Cost-effective screening option
 - Culturally sensitive screening option
 - Easily accessible screening option worldwide




Figure 3: Cervical cancer stages illustration

Product Design Specifications



Figure 4: Women in rural western Africa

- Client Requirements
 - Cost between \$3-5 USD
 - Accessible to demographic
 - Easy to use
 - Provides clear results





Figure. 5 : Incidence and mortality based of cervical cancer on geographic location

- Test should provide at least 70% accurate results
- Biocompatible and non toxic
- Should not be biodegradable
- Product should be single use


Design Matrix (Sample Type)

Designs	#1 Blood		#2 Saliva		#3 Urine	
Categories						
Prior Detection (30)	3/5	18	3/5	18	4/5	24
Ease of Obtaining Usable Sample (25)	4/5	20	3/5	15	4/5	20
Comfort (20)	2/5	8	5/5	20	4/5	16
Ease of Collection (15)	2/5	6	5/5	15	4/5	12
Storage Requirements (10)	2/5	4	5/5	10	4/5	8
Total (100)	56		78		80	


Design 1: Strip Dip



Design 2: Drop Tests



Design 3: Funnel Device



Design Matrix (Collection Device)

Designs	#1 Strip Dip		#2 Drop Test		#3 Funnel Device	
Categories			• —			
Ease of Use (30)	3/5	18	4/5	24	5/5	30
Cost (25)	5/5	25	3/5	15	2/5	10
Ease of Fabrication (20)	5/5	20	4/5	16	2/5	8
Sample Containment (15)	3/5	9	5/5	15	2/5	6
Efficiency (10)	4/5	10	5/5	10	3/5	6
Total (100)	82		80		60	

Josephine Hall

Future Challenges

- Finding an Accessible Biomarker
 - let-7d-3p and miR-30d-5p
 - HPV antibody
 - HPV16/18 E6 oncoprotein
- Cost Considerations
 - Goal of \$3-5 sale price
- Sample collection device modifications






Figure 6: Image of HPV antibody

Figure 7: Image of one step in the process of DNA extraction

Acknowledgements

- Special thanks to:
- -Ms. Kebron Zeygeye
- -Dr. Qian
- -The BME Department

References

- "Human Papillomavirus (HPV)," Planned Parenthood. [Online]. Available: https://www.plannedparenthood.org/learn/stds-hiv-safer-sex/hpv.
 - https://www.plannedparenthood.org/learn/stds-hiv-safer-sex/hpv. [Accessed: 13-Oct-2021].
- K. U. Petry, "HPV and cervical cancer," *Scand. J. Clin. Lab. Invest.*, vol. 74, no. sup244, pp. 59–62, Aug. 2014, doi: 10.3109/00365513.2014.936683. [Accessed: 3-Oct-2021].
- "Cervical cancer," *World Health Organization*. [Online]. Available: https://www.who.int/health-topics/cervical-cancer#tab=tab_1. [Accessed: 23-Sep-2021].
- "Behaviors surrounding oral sex may increase HPV-related cancer risk," *Johns Hopkins Medicine Newsroom*, 14-Jan-2021. [Online]. Available:
 - https://www.hopkinsmedicine.org/news/newsroom/news-releases/b ehaviors-surrounding-oral-sex-may-increase-hpv-related-cancer -risk. [Accessed: 13-Oct-2021].
- "Pap smear," *Mayo Clinic*, 25-Jun-2020. [Online]. Available: https://www.mayoclinic.org/tests-procedures/pap-smear/abo ut/pac-20394841. [Accessed: 16-Sep-2021].
- D. Fields, "HPV and Cervical Cancer," *Medical News*, 27-Feb-2019. [Online]. Available:

https://www.news-medical.net/health/HPV-and-Cervical-Cancer.as px. [Accessed: 13-Oct-2021].

- "New HPV Monoclonal Antibodies," USBIO.NET. [Online]. Available: https://www.usbio.net/promos/hpv-antibodies. [Accessed: 13-Oct-2021].
- Let's Talk Science and A. 30, "DNA extraction," *Let's Talk Science*, 30-Apr-2013. [Online]. Available: https://letstalkscience.ca/educational-resources/backgr ounders/dna-extraction. [Accessed: 13-Oct-2021].
- A. Daponte, G. Michail, A.-I. Daponte, N. Daponte, and G. Valasoulis, "Urine HPV in the context of genital and cervical cancer screening—an update of current literature," *Cancers*, vol. 13, no. 7, p. 1640, 2021.
- R. van den Helder, N. E. van Trommel, A. P. van Splunter, B. I. Lissenberg-Witte, M. C. G. Bleeker, and R. D. M. Steenbergen, "Methylation analysis in urine fractions for optimal CIN3 and cervical cancer detection," *Papillomavirus Research*, vol. 9, p. 100193, 2020
- M. Zheng, L. Hou, Y. Ma, L. Zhou, F. Wang, B. Cheng, W. Wang, B. Lu, P. Liu, W. Lu, and Y. Lu, "Exosomal let-7D-3P and Mir-30d-5p as diagnostic biomarkers for non-invasive screening of cervical cancer and its precursors," *Molecular Cancer*, vol. 18, no. 1, 2019.