

Wearable Simulator For Enhanced Realism

Preliminary Presentation October 15th, 2021





Team Members

Client: Dr. Michael Lohmeier

Section Chief EMS at UW Health, board certified in Emergency Medicine

Advisor: Dr. Krishanu Saha

Team Leader: Francisco O'Neill Rodríguez

Communicator: Sierra Reschke

BSAC: Sami Matloub

BWIG: Jessie Hunt

BPAG: Sami Matloub

Problem Statement

Simulations and mannequins that exist today are able to mock human injuries in impressive and accurate ways, however they are unable to simulate an actual human.



https://www.healthysimulation.com/medical-simulator/

Background

- Purpose and use
- Specific conditions for the vest to simulate
 - Heart conditions: bradycardia, tachycardia, and heart murmurs
 - Lung conditions: Rales, fluid in the lungs, different sounds from right versus left lungs
- Previous design progress

Client Requirements

- No more than \$500
- Must be a wearable simulator of reasonable weight and size to fit the average person
- Can be used 4-5 times a month
- Outputs can be modified during simulation to respond to interventions
- Simulates heart and lung sounds and pulses

Product Design Specifications

- The vest should be comfortable, durable, and easy to use for incoming medical staff
- There will be a single prototype of this medical simulator involved in this project
- Should produce accurate sounds to allow trainees to detect irregularities within heart or lungs
- Vest should be only used within the correct environment
 - Hospitals or medical clinics

Sami Matloub

Existing Devices/Competitors

- → SimMan/3G/Plus
 - Laerdal Medical
- → Simulaids Smart STAT Basic with iPad
 - Universal Medical



- https://laerdal.com/us/products/simulation-training/emergency-care-trauma/simman-3g/
- → HAL S3201 Advanced Multipurpose Patient Simulator
 - Gaumard Scientific



Vest Designs

Condor Ballistic Vest



https://www.amazon.com/Condor-MV-001-Modular-Vest-OliveDrab/dp/B0072 K82TW

- Exterior tough canvas
- Inside padding
- MOLLE design
- Various built-in pockets

Long Sleeve Vest



https://www.wakemakers.com/ronix-supreme-athletic-cut-comp-vest.ht

- Able to run wires up neck and arms
- One piece garment
- Provides trainees with more areas to detect vitals

Add-On Vest




https://uarmprotection.com/product/groin-armor-module/

- Malay system for vest
- Very module and mobile
- Adaptable to wide range of add ons
- Space for extensive wiring

Design Matrix- Vest Design

Design Criteria	Weight	Condor Ballistic Vest		Long Sleeve Vest		Add-On Version	
Durability	30	5/5	30	3/5	18	4/5	24
Output potential	15	4/5	12	4/5	12	5/5	15
Adaptability	15	3/5	9	2/5	6	5/5	15
Feasibility	10	5/5	10	2/5	4	3/5	6
Safety	10	5/5	10	4/5	8	4/5	8
Comfort	5	3/5	3	3/5	3	2/5	2
Aesthetics	5	2/5	2	4/5	4	3/5	3
Sterilizability	5	5/5	5	2/5	2	3/5	3
Cost	5	5/5	5	2/5	3	2/5	2
Total	100		56		42		54

Electronic Designs




Design Matrix- Electronics

Design Criteria	Weight	Dial Controlled Speakers		Wireless/Bluetooth		Wired to Computer	
Ease of Use	25	3/5	15	5/5	25	3/5	15
Output Accuracy	20	4/5	16	4/5	16	4/5	16
Durability	15	3/5	9	4/5	12	3/5	9
Feasability	15	2/5	6	4/5	12	3/5	9
Adaptability	10	3/5	6	3/5	6	3/5	6
Cost	5	3/5	3	3/5	3	4/5	4
Motility	5	4/5	4	5/5	5	2/5	2
Safety	5	3/5	3	4/5	4	2/5	2
Total	100		62		83		63

Future Work

- Visit and tour the ER and Sim Man
- Conduct further research on wireless electronic components
- Order necessary materials
- Begin prototype and testing



References

- US Department of Health and Human Services (1997, February). Good Guidance Practices: Labeling Regulatory Requirements for Medical Devices. HHS Publication FDA89-4203.https://www.fda.gov/files/medical%20devices/published/Labeling---Regulatory-Requirements-for-Medical-Devices-%28FDA-89-4203%29.pdf
- [2] National Institute of Justice, "Current and Future Research on Body Armor," January 1, 2010, https://nij.ojp.gov/topics/articles/current-and-future-research-body-armor
- [3] Bally Ribbon Mills. n.d. Polyethylene | Bally Ribbon Mills. https://www.ballyribbon.com/fibers/ polyethylene
- [4] Total Security Solutions. 2016. Today's Police Body Armor: Overt, Lighter, Smaller. Total Security Solutions https://www.tssbulletproof.com/blog/todays-police-body-armor-overtlighter-smaller
- [5] INACSL Standards Committee (2016, December). INACSL standards of best practice: Simulation SM. Clinical Simulation in Nursing, 12(S), S5-S47. https://www.inacsl.org/INACSL/ document-server/?cfp=INACSL/assets/File/public/standards/SOBPEnglishCombo.pdf.
- [6] Center for Devices and Radiological Health. (n.d.). *Classify your medical device*. U.S. Food and Drug Administration. Retrieved September 23, 2021, from https://www.fda.gov/medical-devices/overview-device-regulation/classify-your-medical-device.
- [7] Simman 3G Plus immersive training for emergency care procedures. Laerdal Medical. (n.d.). Retrieved September 23, 2021, from https://laerdal.com/us/products/simulation-training/emergency-care-trauma/simman-3g/.
- [8] "Simulaids SMART STAT Basic with IPad." Universal Medical, www.universalmedicalinc.com/simulaids-smart-stat-basic-with-ipad.html?campaignid=1049828359.
- [9] Gaumard Scientific. (n.d.). Retrieved October 14, 2021, from https://www.gaumard.com/advanced-hal-s3201.

Acknowledgments

We would like to thank:

- Department of Biomedical Engineering
- Dr. Michael Lohmeier
- Dr. Krishanu Saha