#### Tissue Model of the Epithelial Mesenchymal Trophic Unit



#### DEPARTMENT OF Biomedical Engineering UNIVERSITY OF WISCONSIN-MADISON

Dr. Tracy Jane Puccinelli Dr. Allan Brasier October 6th, 2023

Carley Schwartz (Co-Leader), Elijah Diederich (Co-Leader), Caitriona Treacy (Communicator),

Will Onuscheck (BSAC), Anuraag Shreekanth Belavadi (BWIG), Nick Herbst (BPAG)


### **Problem Statement**

- Chronic lung diseases can cause damage to epithelial tissues of the lungs
  - Pulmonary fibrosis, asthma, and COPD
  - Damage causes the sub-epithelial fibroblasts to increase production
- Currently no scaffolds that accurately model the lung extracellular matrix and its changes due to cell injury
  - Varying mechanical stiffness, porosity, incorporation of collagen and fibronectin, cell adhesive properties
- Dr. Brasier of the UW School of Medicine and Public Health requires such a scaffold that allows for lung epithelial cell culture in an ALI
  - Aim is to study cells in normal and fibrotic ECM conditions



### **Broader Impact**


- Provide a more realistic and accurate 3D cell culture model for lung epithelial cells.
- Improvement to drug Discovery and Development
- Potential for Personalized Medicine
- Application in related fields
- Improve the quality of life for patients with chronic lung diseases.





## Background: Lung Extracellular Matrix

- Lung ECM
  - Collagen, elastin, laminin, and fibronectin [2]
- Function
  - Physical support, cell migration tract
  - Presents and stores growth factors [3]
- Fibroblasts
  - ECM protein production
  - Effector cell for injury repair [4]
- Collagen(I-IV and XVIII)
  - Provides tensile strength, regulates cell adhesion and directs tissue development [5]



#### Figure 2: Lung ECM Diagram [6]



## Background: Cell Culture

Bioprinting

- 3D tissue-like structures with precise spatial control
- Biolnk GelMa-Fibroblasts-Biomolecules

3D cell culture: Hydrogel

- More accurately mimic *in-vivo* ECM
- Allow cell-cell and cell-ECM interactions [7]

Hydrogels:

- Natural
  - Gelatin, Alginate, Collagen
  - Biodegradable, Natural adhesive properties
- Synthetic
  - PEG, PLA
  - Long lasting, replicable [8]






#### Figure 3: Cellink BioX Bioprinter [9]

## **Competing Designs**

Complexity of culture

- 2D Models typically include layers of cells on top of polymer or glass dishes
  - Young's Modulus in 2-4 GPa range 0 [10]
  - Negatively impacts gene expression 0
- **3D Models** 
  - Matrigel 0
    - Derived from Mouse Tumors
    - High variability in batches [11]
  - Human Lung ECM Hydrogels Ο
  - Hyaluronic Acid (HA) Hydrogels Ο
    - Incorporated in PEG hydrogels
    - Free Radical Toxicity

Matrigel



Poorly defined 2D stem-cell culture

Tumour-mimetic

3D cell culture



Chemically defined

- Xenogenic-free Tunable and controllable
- Reproducible
- Broadly applicable





Tissue-mimetic, tunable 3D cell culture



Controlled organoid


assembly

Variable organoid assembly

Figure 4: Competing Scaffold Design [12]



Synthetic scaffolds



cell culture



# Key Design Criteria

- Biochemical Properties
  - Supports cell adhesion and capable of fibroblast encapsulation
  - Porosity
    - Transportation of media through hydrogel
  - Degradation to allow for ECM remodeling
- Mechanical Properties
  - Two types of hydrogels produced: healthy and fibrotic tissue states
    - Healthy Young's Modulus: 2-5 kPa
    - Fibrotic Young's Modulus: ≥16.5 kPa
- Replicability/Ease of fabrication
  - Simple protocol that is easy to follow
  - Accessible materials
    - Pre-characterized Cellink GelMA for easier fabrication



## **Previous Semester**

- Final Design
  - GelMA (50% DOF) hydrogel created by pipetting solution into molds and photocrosslinking with 365 nm wavelength
- Testing
  - Mechanical properties: Rheometry to determine Young's Modulus
    - Lower stiffness reached: 4.08 +/- 0.56 kPa average
    - Higher stiffness reached: 24.2 +/- 9.2 kPa average
  - Cell adhesion: Imaging of hydrogel surface and cell counting
- Lessons Learned
  - Difficult to achieve replicability
    - Variations in mechanical properties between batches following the same protocol
  - GelMA promotes cell viability
    - At lower stiffness there was significantly less cell attachment
      - Incubating in cell media rather than PBS for 24 hours aids in attachment






Figure 5: Image of previous semester cell culture on hydrogel






#### **November Goals**



UNIVERSITY OF WISCONSIN-MADISON



# Budget

- \$5,000 budget given Fall 2022
  - \$679 spent Fall 2022: PEG hydrogel-related materials
  - \$412 spent Spring 2023: GelMA hydrogel-related materials
- Past materials mostly given to client due to project direction changing
- \$3,909 remaining
  - Client provided bioprinter, GelMA bioink, and paraphernalia
  - Potential future purchases
    - Additional GelMA bioink cartridges
    - UV light



#### Acknowledgements

#### Thank you to Dr. Tracy Jane Puccinelli, Dr. John Puccinelli, and Dr. Allan Brasier & team!



#### References

- [1] "Treating chronic obstructive pulmonary disease," McKesson Medical-Surgical,
- https://mms.mckesson.com/content/our-services-solutions/patient-care-engagement/reducing-hospital-readmissions/copd/ (accessed Oct. 4, 2023). [2] E. S. White, "Lung extracellular matrix and fibroblast function," Annals of the American Thoracic Society, Mar-2015. [Online]. Available:
  - https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4430981/
- [3] T. J. Puccinelli, "Biomaterial Scaffold Design," presented in BME 510, University of Wisconsin-Madison, Sep. 14, 2023. [PowerPoint slides].
- [4] G. F; "Fibronectin and wound healing," Journal of cellular biochemistry, 1984. [Online]. Available: https://pubmed.ncbi.nlm.nih.gov/6084665/. [Accessed: 06-Oct-2022].
- [5] "Fibroblasts in lung homeostasis and disease," Thoracic Key, 11-Jan-2017. [Online]. Available: https://thoracickey.com/fibroblasts-in-lung-homeostasis-and-disease].
- [6] P. C. Wijsman, L. H. van Smoorenburg, D. M. de Bruin, J. T. Annema, H. A. M. Kerstjens, O. M. Mets, M. van den Berge, P. I. Bonta, and J. K. Burgess, "Imaging the pulmonary extracellular matrix," Current Opinion in Physiology, vol. 22, p. 100444, Aug. 2021
- [7] H. Chamkouri and M. Chamkouri, "A Review of Hydrogels, Their Properties and Applications in Medicine," Biomedgrid, 03-Feb-2021. [Online]. Available: https://biomedgrid.com/fulltext/volume11/a-review-of-hydrogels-their-properties-and-applications-in-medicine.001682.php
- [8] "2D vs 3D Creative Biolabs," Cell Culture: 2D Vs 3D Creative Biolabs. [Online]. Available: https://www.creative-biolabs.com/adc/2d-vs-3d.htm
- [9] "Bio XTM 3D bioprinter," CELLINK, https://www.cellink.com/bioprinting/bio-x-3d-bioprinter/ (accessed Oct. 4, 2023).
- [10] D. T. Butcher, T. Alliston, and V. M. Weaver, "A tense situation: Forcing tumour progression," Nature reviews. Cancer, 01-Feb-2009. [Online]. Available: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2649117/. [Accessed: 22-Feb-2023].
- [11] S. Kim, S. Min, Y. S. Choi, S.-H. Jo, J. H. Jung, K. Han, J. Kim, S. An, Y. W. Ji, Y.-G. Kim, and S.-W. Cho, "Tissue extracellular matrix hydrogels as alternatives to Matrigel for culturing gastrointestinal organoids," Nature News, 30-Mar-2022. [Online]. Available: https://www.nature.com/articles/s41467-022-29279-4. [Accessed: 23-Feb-2023].
- [12] E. A. Aisenbrey and W. L. Murphy, "Synthetic alternatives to Matrigel," Nat Rev Mater, vol. 5, no. 7, Art. no. 7, Jul. 2020, doi: 10.1038/s41578-020-0199-8.

