

Far-UVC Light in Clinical Settings

Client: Dr. Ernesto Brauer

Advisor: Randolph Ashton

Team Members: Parker Esswein, Lars Krugel, Tyler Linderman,
Derick Peters, Vanessa Obrycki and Draeson Marcoux





Our Client

- Dr. Ernesto Brauer
- Critical Care Physician at Aurora St. Luke's Medical Center
- Personal use of UVC disinfection device
- Very hopeful for future applications of far-UVC in clinical settings




<https://www.doximity.com/pub/ernesto-brauer-md>



<https://www.grainger.com/product/AIR-SCIENC-E-Stainless-Steel-UV-Box-Sterilization-18AX46>

Problem Statement and Requirements


- Design and fabricate a Far-UVC light device to be used in clinical settings
 - Full room coverage
 - 99.9% effective against common household bacteria (fecal patina)
- Safe for constant human exposure
- Determine exposure time, distance and intensity relationship






Background Research

- 222 nm vs. 254 nm wavelengths (Buonanno, 2020)
 - Penetration - Living cells in the skin or eyes
- Two current studies
 - Columbia Medical Center Research (Buonanno, 2020)
 - Human Lung Cells
 - HCoV-229E and HCoV-OC43
 - 22 cm = 99.9% effective in 25 min at 1.7 mJ/cm²
 - Hiroshima University Research Group (Hiroshima, 2020)
 - Solution of SARS-CoV-2 on a plate
 - 24 cm = 99.7% effective in 30 seconds at 3.6mJ/cm²
- Both in vitro
- Minimal distances



<https://marktchopto.com/technical-articles/understanding-ultraviolet-led-applications-and-precautions/>



Ushio - 222-nm KrCl excimer lamp

<https://www.ushio.com/product/care222-filtered-far-uv-c-excimer-lamp-module/>


Project Design Specifications

- Function:
 - Kill bacteria efficiently
 - Use of Far-UVC wavelength
- Requirements:
 - Performance
 - Safety
 - Life/Durability
 - Environment
 - Physical Characteristics
- Production:
 - Cost
 - Quantity
- Other:
 - Standards
 - Competition
 - Patient/Client Preferences




FULED Overhead Light



FULED



2-in-1 Air Purifier and Far-UVC Sanitation



Mobile Light Cart - easy access





Design Matrix Criteria

1. Efficacy
2. Coverage
3. Safety
4. Ease of Fabrication
5. Cost
6. Energy Usage
7. Durability
8. Ease of Use



<https://hydrobuilder.com/learn/grow-light-coverage-calculator/>


Far-UVC Device Designs

Design Criteria (Weight)	Mobile Light Cart- Easy Access		FULED Overhead Light		2-in-1 Air Purifier and Far-UVC Sanitation	
Efficacy*(25)	4/5	20	5/5	25	4/5	20
Coverage/Square Footage (20)	3/5	12	5/5	20	2/5	8
Safety (15)	4/5	12	2/5	6	4/5	12
Ease of Fabrication (15)	3/5	9	2/5	6	4/5	12
Cost(10)	3/5	6	2/5	4	2/5	4
Energy Usage (5)	4/5	4	2/5	2	2/5	2
Durability (5)	4/5	3	5/5	5	2/5	2
Ease of Use (5)	3/5	3	5/5	5	4/5	4
Total (100)		66		73		58

Design Matrix

Future Work

- Begin developing the optimal design
- Use calculations to determine:
 - Efficacy
 - Distance light travels and coverage
 - Effectiveness on different materials
 - Various shapes, angles and objects in the room
 - Safety
 - Intensity
 - Continued exposure
- Test calculations and design(s)





Acknowledgements

We would like to thank Dr. Ashton, Dr. Rogers, and our client, Dr. Brauer, for being outstanding mentors and support systems for us throughout the development of this project.

References

Brenner, David. “Can Light Stop the Coronavirus?” *TED*, www.ted.com/talks/david_brenner_can_light_stop_the_coronavirus.

Buonanno, M., Welch, D., Shuryak, I. *et al.* Far-UVC light (222 nm) efficiently and safely inactivates airborne human coronaviruses. *Sci Rep* 10, 10285 (2020). <https://doi.org/10.1038/s41598-020-67211-2>

“Study Shows First Proof That a Safer UV Light Effectively Kills Virus Causing COVID-19.” 広島大学, www.hiroshima-u.ac.jp/en/news/60119.