

Client: Dr. Frederick Kelcz


Advisors: Dr. Lonie Salkowski and Dr. Beth Meyerand

Gopika Senthilkumar
Kevin Fantl
Colin Schrof
Alex Henry
Kari Borowski

# **Background - Current Method**














#### **Problem Statement**





Standardize localization technique



- Clinically effective
- Accurate Perpendicular approach
- Easy and quick to use
- Safe for patients
- Safe for physicians
- Removable/ Radiotranslucent

#### **Justification**



- Standardization Quality Control and Save Cost
- Easier to learn procedure
- More comfortable experience for patients
  - Time
  - Reduced radiation
  - Not freehand
- Efficiency
- Technology may be applied in other situations


#### **Last Semester**



- Finalized design of guide
- Designed and validated new Breast Model
- Perfected and practiced study method
- Started IRB approved clinical testing three subjects
- Started drafting the publication

# **Final Design of Guide**






#### **New Breast Model Validation**







# **Study Design: Overview**

Step 1:

Standard localization procedure on model

Step 2:

Training with needle guide design (30 mins)

Step 3:

Localization using needle on model

Step 4:

Post-procedure survey



### **Study Design: Participants**

- Breast Imaging Fellows
- Physicians
- Residents on Clinical Service
- Any staff authorized to perform localizations











- 1. Validate that our device can accurately guide the needle to its intended location
- 2. More efficient and standardized than the traditional technique
- 3. User friendly design. Does it feel comfortable to use?
- 4. Better for new learners? Other Correlations?

#### **Final Design Testing Results**









Top-image Orthogonal-view

Free-Hand Localization







Top-image Orthogonal-view Localization with Guide

#### **Clinical Testing Results**








- Reduced localization time for ¾ subjects
- Maintained accuracy \*n=3

#### **Draft of Publication**


- Goal: Introduce a validated, Novel Needle
   Guide for standardizing and improving
   the efficiency of Preoperative Wire-Guided
   Localizations of Breast Lesions
  - Maybe:
    - Easier for new learners
    - correlation with current skill level
- Potential Journal: IEEE Journal of Translational Engineering in Health and Medicine
  - Maybe higher impact based on results


# **Semester Goals**



#### **Data Analysis Plan**

- Difference in time from initial needle puncture to final removal of needle
- 2. Perpendicular puncture
- Distance between clip and needle
- 4. Number of corrections during procedure
- Other correlations from Survey
- 6. Statistical Analysis









- Procedure with device is faster while maintaining/improving accuracy
- 2. Gets perpendicular puncture every time
- 3. Device is comfortable to use

\*Other client requirements already met





- 1. Provisional Patent
- 2. Find faster, cheaper, more sterile manufacturing method.
- 3. Each guide packaged in sterile bags one per patient use.
- 4. Include pictographic user guide maybe send representatives to train users
- 5. Include safety information for patients and physicians.
- 6. Approach Hologic with our device.



- Fall Semester Budget: \$200.00
- Spring Semester Budget: <\$100.00

Description	Part #	Source	Quantity	Cost (Per Unit)
Humimic Phantom gel	852844007819	Humimic Medical	1	\$46.05
4x4x4 tin mold for phantom gel	111-4682432-464	Amazon.com	1	\$29.99
Humimic Phantom gel	2955	Humimic Medical	1	\$35.11
Flesh tone, opague pigment dye	111-5042561-625	Amazon.com	1	\$11.50
Printed Design presentation poster	N/A	UW-Madison	1	\$48.00
		9		
	7. P			
	8			200
	*			·

#### **Acknowledgements**

Dr. Lonie Salkowski

Dr. Beth Meyerand

Dr. Frederick Kelcz

BME Design Faculty

**UW Makerspace Faculty** 

