

Tissue Model of the Epithelial Mesenchymal Trophic Unit



DEPARTMENT OF
Biomedical Engineering
UNIVERSITY OF WISCONSIN-MADISON

Prof. Tracy Jane Puccinelli

Dr. Allan Brasier

February 9th, 2024

Carley Schwartz (Co-Leader), Elijah Diederich (Co-Leader), Caitriona Treacy (Communicator),

Will Onuscheck (BSAC), Anuraag Shreekanth Belavadi (BWIG), Nick Herbst (BPAG)

Problem Statement

- Chronic lung diseases can cause damage to epithelial tissues of the lungs
 - Pulmonary fibrosis, asthma, and COPD
 - Damage causes the sub-epithelial fibroblasts to increase production
- Currently no scaffolds that accurately model the lung extracellular matrix and its changes due to cell injury
 - Varying mechanical stiffness, porosity, incorporation of collagen and fibronectin, cell adhesive properties
- Our client, Dr. Brasier of the UW School of Medicine and Public Health, requires such a scaffold that allows for lung epithelial cell culture in an ALI
 - Aim is to study cells in normal and fibrotic ECM conditions



Broader Impact

- COPD is 3rd leading cause of death in the world [1]
 - In the US, over 120,000 people die yearly of COPD
- A tissue model would...
 - Contribute to more thorough understanding of chronic lung diseases
 - Aid in the development of individualized treatments/therapies



Figure 1: Graphical depiction of chronic lung diseases [1]



Background: Lung Extracellular Matrix

- Lung ECM
 - Collagen, elastin, laminin, and fibronectin [3]
- Function
 - Physical support, cell migration tract
 - Presents and stores growth factors [4]
- Fibroblasts
 - ECM protein production
 - Effector cell for injury repair [5]
- Collagen(I-IV and XVIII)
 - Provides tensile strength, regulates cell adhesion and directs tissue development [6]




Figure 2: Lung ECM Diagram [7]



Background: Cell Culture

Bioprinting

- 3D tissue-like structures with precise spatial control
- BioInk - GelMa-Fibroblasts-Biomolecules

3D cell culture: Hydrogel

- More accurately mimic *in-vivo* ECM
- Allow cell-cell and cell-ECM interactions [8]

Hydrogels:

- Natural
 - Gelatin, Alginate, Collagen
 - Biodegradable, Natural adhesive properties
- Synthetic
 - PEG, PLA
 - Long lasting, replicable [9]



Figure 3: Cellink BioX Bioprinter [10]



Competing Designs

- 2D Models typically include layers of cells on top of polymer or glass dishes
 - Young's Modulus in 2-4 GPa range [11]
 - Negatively impacts gene expression
- 3D Models
 - Matrigel
 - Derived from Mouse Tumors
 - High variability in batches [12]
 - Human Lung ECM Hydrogels
 - Hyaluronic Acid (HA) Hydrogels
 - Incorporated in PEG hydrogels
 - Free Radical Toxicity




Figure 4: Competing Scaffold Design [12]



Key Design Criteria

- Biochemical Properties
 - Supports cell adhesion and capable of fibroblast encapsulation
 - Porosity
 - Transportation of media through hydrogel
 - Degradation to allow for ECM remodeling
- Mechanical Properties
 - Two types of hydrogels produced: healthy and fibrotic tissue states
 - Healthy Young's Modulus: 2-5 kPa
 - Fibrotic Young's Modulus: ≥ 16.5 kPa
- Replicability/Ease of fabrication
 - Simple protocol that is easy to follow
 - Accessible materials
 - Pre-characterized Cellink GelMA for easier fabrication



Previous Semester

- Final Design
 - Pipette-based GelMA hydrogel for testing and comparison
 - Bioprinted GelMA hydrogel as final design
- Testing
 - Rheometry to determine Young's Modulus of pipette-based hydrogels
 - Varied setting time at 4 °C and UV-crosslinking time
 - LIVE/DEAD staining and gel imaging for preliminary cell viability testing
 - Staining results were unusable, need to identify new method of quantifying cell viability
 - Bioprinter Optimization
 - Changed equilibrium time, heating and printing temperatures, extrusion pressure, infill density, and UV time

Pipette-Based Hydrogels	
Condition (4 °C, UV)	Stiffness (kPa)
3 min, 5 min	4.2 ± 0.92
5 min, 5 min	6.63 ± 2.6
10 min, 5 min	13.24 ± 2.8
5 min, 7 min	49.65 ± 22
5 min, 10 min	277 ± 155

Figure 5: Mechanical testing results for pipette-based hydrogel batches over course of semester..



Previous Semester

- Lessons Learned
 - Difficult to achieve replicability
 - Variations in mechanical properties between batches following the same conditions
 - Keep a careful eye on the shelf life of the GelMA bioink
 - Write up protocols before conducting experiments
 - Have a statistically significant number of samples ($n \geq 4$)
 - Have literature-based evidence to back up the protocols
 - Bioprinter parameters need to be adjusted more
 - UV time, extrusion pressure, and infill density especially



Semester Timeline

Task	Jan	Feb					March					April				May		
	26	2	9	16	23	28	1	8	15	22	5	12	19	26	1	3	10	
Project Goals																		
Pipette-Based Hydrogel Characterization																		
Bioprinted Hydrogel Protocol																		
Consistent Bioprinted Hydrogel Fabrication																		
Bioprinted Hydrogel Characterization																		
Fibroblast Encapsulation																		
Fibroblast Viability Testing																		
Epithelial Cell Culture w/ ALI																		

Task	Jan	Feb					March				April				May		
	26	2	9	16	23	28	1	8	15	22	5	12	19	26	1	3	10
Deliverables																	
Progress Reports																	
Journal Selection																	
Preliminary Presentation																	
Preliminary Report																	
Preliminary Notebook																	
Show and Tell																	
Executive Summary Draft																	
Executive Summary																	
Final Poster																	
Final Report																	
Final Notebook																	
Client Evaluation																	



Budget

- \$5,000 budget given Fall 2022
 - \$679 spent Fall 2022: PEG hydrogel-related materials
 - \$412 spent Spring 2023: GelMA hydrogel-related materials
 - \$395 spent Fall 2023: Bioprinting-related materials
- Past materials mostly given to client due to project direction changing last semester
- \$3,514 remaining
 - Client provided bioprinter and paraphernalia
 - Future purchases will primarily be more bioink cartridges and bioprinter paraphernalia
 - \$325 for 3 x 3mL cartridges
 - \$990 for 10 x 3mL cartridges
 - \$85 for a 50-count pack of 22 gauge nozzles



References

- [1] "Chronic obstructive pulmonary disease (COPD)." Accessed: Dec. 06, 2023. [Online].
- [2] E. S. White, "Lung extracellular matrix and fibroblast function," *Annals of the American Thoracic Society*, Mar-2015. [Online]. Available: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4430981/>
- [3] T. J. Puccinelli, "Biomaterial Scaffold Design," presented in BME 510, University of Wisconsin-Madison, Sep. 14, 2023. [PowerPoint slides].
- [4] G. F. "Fibronectin and wound healing," *Journal of cellular biochemistry*, 1984. [Online]. Available: <https://pubmed.ncbi.nlm.nih.gov/6084665/>. [Accessed: 06-Oct-2022].
- [5] "Fibroblasts in lung homeostasis and disease," *Thoracic Key*, 11-Jan-2017. [Online]. Available: <https://thoracickey.com/fibroblasts-in-lung-homeostasis-and-disease>.
- [6] P. C. Wijsman, L. H. van Smoorenburg, D. M. de Bruin, J. T. Annema, H. A. M. Kerstjens, O. M. Mets, M. van den Berge, P. I. Bonta, and J. K. Burgess, "Imaging the pulmonary extracellular matrix," *Current Opinion in Physiology*, vol. 22, p. 100444, Aug. 2021
- [7] H. Chamkouri and M. Chamkouri, "A Review of Hydrogels, Their Properties and Applications in Medicine," *Biomedgrid*, 03-Feb-2021. [Online]. Available: <https://biomedgrid.com/fulltext/volume11/a-review-of-hydrogels-their-properties-and-applications-in-medicine.001682.php>
- [8] "2D vs 3D - Creative Biolabs," *Cell Culture: 2D Vs 3D - Creative Biolabs*. [Online]. Available: <https://www.creative-biolabs.com/adc/2d-vs-3d.htm>
- [9] "Bio XTM - 3D bioprinter," *CELLINK*, <https://www.cellink.com/bioprinting/bio-x-3d-bioprinter/> (accessed Oct. 4, 2023).
- [10] D. T. Butcher, T. Alliston, and V. M. Weaver, "A tense situation: Forcing tumour progression," *Nature reviews. Cancer*, 01-Feb-2009. [Online]. Available: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2649117/>. [Accessed: 22-Feb-2023].
- [11] S. Kim, S. Min, Y. S. Choi, S.-H. Jo, J. H. Jung, K. Han, J. Kim, S. An, Y. W. Ji, Y.-G. Kim, and S.-W. Cho, "Tissue extracellular matrix hydrogels as alternatives to Matrigel for culturing gastrointestinal organoids," *Nature News*, 30-Mar-2022. [Online]. Available: <https://www.nature.com/articles/s41467-022-29279-4>. [Accessed: 23-Feb-2023].
- [12] E. A. Aisenbrey and W. L. Murphy, "Synthetic alternatives to Matrigel," *Nat Rev Mater*, vol. 5, no. 7, Art. no. 7, Jul. 2020, doi: 10.1038/s41578-020-0199-8.

