

Vaginal Self-Swab to Minimize Contact Contamination

Sara Morehouse (Leader), Cherry Qiu (Communicator), Katherine Kafkis (BWIG and BSAC), and Adam Berdusco (BPAG) BME 301 Advisor: Dr. Megan McClean TA: Stephen Foley Client: Dr. Jean Riquelme 03/01/24

Outline

- Problem Statement
- Background Information
- Product Design Specifications
- Specific Semester Goals
- Design 1 Modified Plunger
- Design 2 Snap On
- Design 3 Pull Back
- Design Matrix
- Future Work
- Acknowledgements & Questions
- References



Problem Statement

Current STI self-testing methods contaminate the testing environment

- 84% prefer self-swabs to gynecological exams [1]
- 67% of women received a false positive result due to contamination [2]

Contamination Methods

- Insertion into and swabbing of vaginal canal
- Transfer to the media tube



Background Information

Client: Dr. Jean Riquelme

- Family medicine physician based in Madison
- Requesting a vaginal self-swab device for use in Chlamydia screening
- Promote universal STI testing



Sara

Background Information



Fig. 2: Aptima Multitest swab instructional image [5]

- STIs are under-tested especially in young women
 - 1 in 4 sexually active adolescent females in US have an STI
 [4]
 - Asymptomatic
 - Barriers to testing [4]
 - Long term complications if untreated [1]
- Nucleic Acids Amplification Test [1]
- Existing designs
 - Aptima test by Hologic [5]
 - Patient RNA remained on testing room surfaces [2]
 - Require more thorough cleaning



Product Design Specifications

<u>Goal:</u> Create a vaginal self-swab device that reduces contamination of the testing room while prioritizing patient comfort in order to promote universal testing.

- User-friendly
- Housing/deployment and sealing mechanisms
 - Release of transport media [1]
- Head of swab must insert 5 cm into the vagina [3]
- Biocompatible and non-toxic materials
 - Non-cotton fiber swab (Dacron) [7]
 - Universal transport media
 - Biocompatible plastics: PLA, PP, etc.
- Overall length under ~17 cm
- Manufactured with 3D-printing
- Budget: \$500

Fig. 3 Materials provided with current self-swab tests [5]

Katherine

Specific Semester Goals

Leaking - Transport Media

- Class 1 medical device [8]
 - Required for sample preservation [1]
- Must be contained within the device

Sealing

• Induction-sealed thin film within media container (cap)

Ergonomics/Aesthetics

- Thin in diameter, short in length
- Threading on inside of device
- Arrows or colors incorporated for further direction

Sustainability - Single Use

• Compostable plastics: PLA, PBS, PBAT mixtures [9]



Fig.4 Previous device with areas of improvement indicated.

Katherine

Design Idea 1 - Modified Plunger Design



• Similar to a tampon

- Plunger holds swab and provides 6 cm of vertical motion
- Body contains the swab when it is not inserted into the vagina
- Two o-rings prevent leakage of media
- Cap has location for thin film attachment
- Threads on the interior of device for aesthetics

Adam

Design Idea 2 - Snap On Design



Fig.6 Full and sectional views of the snap on design

- Intuitive design that allows for swabbing with one hand
- A second locking mechanism allows the swab to be submerged in media
- Media is housed in the bottom of the base
- The device will not leak if kept upright

Adam

Design Idea 3 - Pull Back Design



- Twisting lock mechanism allows for easy access to the swab
- The top sheath protects the swab from contamination
- Media is housed in the top of the cap (not shown)
- Room for two o-rings will help prevent leaking

Adam

Design Matrix

Weight	Modified Plunger		Snap On		Pull Back	
30	5/5	30	3/5	18	5/5	30
25	4/5	20	5/5	25	4/5	20
15	4/5	12	5/5	15	3/5	9
10	1/5	2	5/5	10	3/5	6
10	5/5	10	4/5	8	4/5	8
5	5/5	5	5/5	5	5/5	5
5	5/5	5	5/5	5	5/5	5
100	84		86		83	
	Weight 30 25 15 10 5 5 100 30 30 30 30 30 30 3	Weight Modified Plus 30 5/5 30 5/5 25 4/5 15 4/5 10 1/5 5 5/5 5 5/5 100 84	Weight Modified Plunger 30 5/5 30 30 5/5 30 25 4/5 20 15 4/5 12 10 1/5 2 5 5/5 5 5 5/5 5 100 84	Weight Modified Plunger Snap Or 30 5/5 30 3/5 30 5/5 30 3/5 25 4/5 20 5/5 15 4/5 12 5/5 10 1/5 2 5/5 5 5/5 5 5/5 5 5/5 5 5/5 100 84 86	Weight Modified Plunger Snap On 1	Weight Modified Plunger Snap On Pull E 30 $5/5$ 30 $3/5$ 18 $5/5$ 30 $5/5$ 30 $3/5$ 18 $5/5$ 25 $4/5$ 20 $5/5$ 25 $4/5$ 15 $4/5$ 12 $5/5$ 15 $3/5$ 10 $1/5$ 2 $5/5$ 10 $3/5$ 10 $5/5$ 10 $4/5$ 8 $4/5$ 5 $5/5$ 5 $5/5$ $5/5$ $5/5$ 5 $5/5$ 5 $5/5$ $5/5$ $5/5$ 100 84 86 86

Future Work

- Redesign prototypes due to limitations of testing equipment
- 3D print designs for fabrication.



Fig.9 Automated diagnostic device that will be used to test the swabs

Cherry

Thank you!

Special thanks to Dr. Riquelme, Dr. McClean, and Stephen Foley for their guidance and input!

Questions?



References

[1]M. Muljadi, C.-M. Cheng, C.-Y. Yang, T.-C. Chang, and C.-J. Shen, "A pilot clinical validation study of a self-collected vaginal swab device for the detection of chlamydia trachomatis in women," *Frontiers in Bioengineering and Biotechnology*, vol. 10, Oct. 2022, doi: https://doi.org/10.3389/fbioe.2022.1008761.

[2]M. Toepfe, B. Hermann, M. Sansone, C. Lilja, and P. Nolskog, "Environmental contamination by Chlamydia trachomatis RNA can cause false-positive test results in clinical samples," *Sexually Transmitted Diseases*, vol. Publish Ahead of Print, Oct. 2020, doi: https://doi.org/10.1097/olq.00000000001323.

[3]"Self-Collected Vaginal Swabs for Gonorrhea and Chlamydia." Accessed: May 05, 2020. [Online]. Available: https://epi.dph.ncdhhs.gov/cd/lhds/manuals/std/labtesting/selfcollectedswabs.pdf

[4]R. J. Steiner, S. L. Michael, J. E. Hall, L. C. Barrios, and L. Robin, "Youth Violence and Connectedness in Adolescence: What Are the Implications for Later Sexually Transmitted Infections?," *Journal of Adolescent Health*, vol. 54, no. 3, pp. 312-318.e1, Mar. 2014, doi: https://doi.org/10.1016/j.jadohealth.2013.09.008.

[5]"Aptima® Multitest Swab | Hologic," *www.hologic.com*. https://www.hologic.com/hologic-products/collection-devices/aptima-multitest-swab (accessed Sep. 12, 2023).

[6]M. J. Gibert, C. Sánchez-Contador, and G. Artigues, "Validity and acceptance of self vs conventional sampling for the analysis of human papillomavirus and Pap smear," *Scientific Reports*, vol. 13, no. 1, p. 2809, Feb. 2023, doi: https://doi.org/10.1038/s41598-023-29255-y.

[7]R. N. Kashapov and A. N. Tsibin, "Comparison of the Physical Properties and Effectiveness of Medical Swabs for Sampling Biomaterials," *Biomedical Engineering*, vol. 55, no. 4, pp. 289–293, Nov. 2021, doi: <u>https://doi.org/10.1007/s10527-021-10120-z</u>.

[8] Transport culture medium, 21 C.F.R. § 866.2390 (2023).

[9] G. Cappiello, C. Aversa, and M. Barletta, "Design of compostable materials for the manufacturing of flexible tampon applicators," Procedia CIRP, vol. 110, pp. 342–347, 2022, doi: https://doi.org/10.1016/j.procir.2022.06.061.