



Inconspicuous Ankle Foot Orthosis (AFO) for Teen

Team Rise and Stride

Lucy Hockerman (Team Leader)

Presley Hansen (BSAC)

Maddie Michels (Communicator)

Kate Hiller (BPAG)

Sadie Rowe (BWIG)

Client **Debbie Eggleston**

Advisor **Dr. John Puccinelli**

TA **Lizzie Maly**

February 21, 2025

Overview

1. Background
2. Problem Statement
3. Previous Work
4. Competing Designs
5. Product Design Specifications
6. Rigid Support Design Matrix
7. Materials Design Matrix
8. Final Design
9. Future Work & Testing
10. References & Acknowledgements



Background

Facioscapulohumeral muscular dystrophy (FSHD):

A rare neuromuscular condition stemming from a genetic mutation that causes progressive muscle weakness [1]

Client background: Debbie Eggleston, physical therapist and activist

Patient: High school student with FSHD

Global Impact: Three to five people out of every 100,000 have FSHD [1]




Figure 1: Ankle Inversion Diagram [2]




Figure 2: Ankle Dorsiflexion Diagram [3]

Problem Statement

Project Description: This project aims to raise awareness for FSHD and explore the benefits of discrete Ankle Foot Orthoses (AFOs) for individuals with this progressive muscle weakness. The team aims to design a comfortable, flexible, and discrete brace tailored for a young individual diagnosed with FSHD to facilitate natural gait.



Previous Work

- Compression sock with gel pads
- Foot brace with CF-PLA piece
- Adjustable bungee cord system
- Gait Analysis
 - Pros: comfortable, no negative impact on gait, and aids dorsiflexion
 - Cons: not easily removable and rigid support



Figure 3: Final Design on Mold of Patient



Figure 4: Final Design on Team Member's Foot for Testing

Competing Designs



Figure 5: Supra Malleolar Orthosis (SMO) [4]



Figure 6: Jointed AFO [5]



Figure 7: Passive Dynamic-AFO [6]



Figure 8: Variable Stiffness Orthosis [7]

Product Design Specifications

- Patient-specific dimensions and comfort needs
- Mimic normal gait and allow more than 30° of motion from neutral ankle position [8]
- Support dorsiflexion, provide resistance 5-10 Nm per 10° of plantarflexion [8]
- Prevent inversion angles greater than 25° [9]
- Rigid support resists torsional forces ± 30 Nm and withstand 266 N of force [10]
- Discrete design for day-to-day use, including horseback riding




Figure 9: Foot Drop With Negative Pitch Angle [11]



Figure 10: Ankle Inversion vs. No Inversion [12]

Rigid Support Designs

Pivot Pro




Figure 11: Pivot Pro Design
Medial and Posterior Views

Calf Hugger




Figure 12: Calf Hugger Front
and Lateral Views

We Support U







Figure 13: We Support U
Medial and Posterior Views

Criteria	 Design 1: Pivot Pro		 Design 2: Calf Hug		 Design 3: We Support U	
	Raw Score	Weighted Score	Raw Score	Weighted Score	Raw Score	Weighted Score
Dorsiflexion range of motion (20)	3/5	12/20	2/5	8/20	5/5	20/20
Mediolateral support (20)	3/5	12/20	3/5	12/20	4/5	16/20
Ease of user assembly (15)	4/5	12/15	4/5	12/15	3/5	9/15
Comfort (15)	2/5	6/15	3/5	9/15	4/5	12/15
Discreteness (10)	3/5	6/10	4/5	8/10	2/5	4/10
Ease of Fabrication (10)	1/5	2/10	2/5	4/10	4/5	8/10
Cost (5)	3/5	3/5	4/5	4/5	3/5	3/5
Safety (5)	4/5	4/5	5/5	5/5	5/5	5/5
Total	57/100		62/100		77/100	

Rigid Support Matrix

Table 1: Rigid Support Design Matrix

Material Considerations

Carbon Fiber Reinforced PLA composite (CF-PLA)



Figure 14: CF-PLA material [13]

- Used in current design
- 3D printed, compatible with 3D scanning
- 470 MPa flexural strength [14]

Fiberglass Plaster



Figure 15: CFN Medical Fiberglass Casting Tape [15]

- Multiple color options, easily accessible
- Simple water activation application process [16]
- 50 MPa flexural strength when wrapped at least twice [17]

Thermoplastics




Figure 16: Thermoplastic sheet vacuum-formed onto objects [18]

- Most common method for AFO fabrication [19]
- Heated thermoplastic sheets vacuum-formed onto mold [20]
- Flexural strength [21]:
 - Polypropylene (PP): 10-20 MPa
 - Polyethylene (HDPE): 10-50 MPa

Criteria	Carbon Fiber reinforced PLA composite (CF-PLA)		Fiberglass Plaster		Thermoplastics	
	Raw Score	Weighted Score	Raw Score	Weighted Score	Raw Score	Weighted Score
Strength/rigidity (30)	5/5	30/30	4/5	24/30	4/5	24/30
Ease of Fabrication (20)	4/5	16/20	5/5	20/20	1/5	4/20
Cost (20)	5/5	20/20	3/5	12/20	4/5	16/20
Safety (20)	5/5	20/20	3/5	12/20	5/5	20/20
Environmental Impacts (10)	5/5	10/10	4/5	8/10	2/5	4/10
Total	96/100		76/100		68/100	

Materials Design Matrix

Table 2: Materials Design Matrix

Final Design




Figure 17: Rigid Support Final Design (We Support U)

Rigid Support

- Free range of dorsiflexion
- Supportive
- Simple assembly

Material

- High strength
- Inexpensive & Accessible
- Environmentally conscious



Figure 18: Carbon Fiber Reinforced PLA Composite (CF-PLA) [13]

Future Work

Prototyping

- Model client's foot
 - 3D scanning cast
 - Silicone molding
- Post processing for 3D printing
 - Alternative: fiberglass plaster



Figure 19:
Provided cast
of client's foot

Testing & Compatibility

- Team & Client evaluation
 - Comfort
- IMU and MoCap Testing via OptiTrack
 - Quantify angle inversion
 - Inversion $< 25^\circ$
 - Track Dorsiflexion
 - ROM $> 30^\circ$ from resting
- Finite Element Analysis (FEA)
 - Force: 266 N
 - Torsion: ± 30 Nm

Acknowledgements

Advisor

Dr. John Puccinelli



TA

Lizzie Maly




Client

Debbie Eggleston



Design Engineer

Jesse Darley



References

- [1] Facioscapulohumeral muscular dystrophy (Fshd). (n.d.). Muscular Dystrophy UK. Retrieved October 3, 2024, from <https://www.muscular dystrophyuk.org/conditions/a-z/facioscapulohumeral-muscular-dystrophy-fshd/>
- [2] "Ankle Joint Eversion and Inversion," *KinX Learning*. <https://kinxlearning.com/blogs/news/ankle-joint-eversion-and-inversion>
- [3] B. Guide, "Best Ankle Sleeves," *Buyersguide.org*, 2025. <https://buyersguide.org/ankle-sleeve/t/best?Country=US&msclkid=73703d6fbb5814edf29b5ad983480b7e&m=e&d=c&c=74148629506956&oid=kwd-74148644213508:loc-190&qs=ankle%20compression%20sleeve&lp=105696&li=8nw=o&nfs=1&tddid=10864757> (accessed Feb. 19, 2025).
- [4] "surestep smo," Surestep. <https://surestep.net/products/surestep-smo/>
- [5] "Professionals," Foot Performance Center. Accessed: Feb. 19, 2025. [Online]. Available: <http://footperformance.com/professional/>
- [6] M. C. Faustini, R. R. Neptune, R. H. Crawford, and S. J. Stanhope, "Manufacture of Passive Dynamic Ankle-Foot Orthoses Using Selective Laser Sintering," *IEEE Trans. Biomed. Eng.*, vol. 55, no. 2, pp. 784–790, Feb. 2008, doi: 10.1109/TBME.2007.912638.
- [7] L. van Noort, N. Van Crey, E. J. Rouse, I. Martinez-Caballero, E. H. F. van Asseldonk, and C. Bayón, "A usability study on the inGAIT-VSO: effects of a variable-stiffness ankle-foot orthosis on the walking performance of children with cerebral palsy," *Journal of NeuroEngineering and Rehabilitation*, vol. 21, no. 1, p. 132, Aug. 2024, doi: 10.1186/s12984-024-01433-7.
- [8] M. Alam, I. A. Choudhury, and A. B. Mamat, "Mechanism and Design Analysis of Articulated Ankle Foot Orthoses for Drop-Foot," *Sci. World J.*, vol. 2014, no. 1, p. 867869, 2014, doi: 10.1155/2014/867869.
- [9] "Ankle Foot Orthosis (AFO) stiffness design for mitigation of ankle inversion injury | IEEE Conference Publication | IEEE Xplore." Accessed: Feb. 02, 2025. [Online]. Available: <https://ieeexplore.ieee.org/abstract/document/8692490>
- [10] "Watch Your Step: Understanding the Importance of Your Feet in the Kinetic Chain – Foundation for Chiropractic Progress," Foundation for Chiropractic Progress - Educating the public to the benefits of chiropractic care. Accessed: Feb. 02, 2025. [Online]. Available: <https://www.f4cp.org/watch-your-step-understanding-the-importance-of-your-feet-in-the-kinetic-chain/>
- [11] Foot, "Polio Network," Polio Network, Feb. 2023, doi: <https://doi.org/101387137/Foot%2BDrop>.
- [12] "Access Library Resource — UW—Madison Libraries," Wisc.edu, 2025. <https://www-sciencedirect-com.ezproxy.library.wisc.edu/science/article/pii/S1350453316301230?via%3Dihub> (accessed Feb. 18, 2025).
- [13] "Amazon.com: FLASHFORGE Carbon Fiber PETG Filament 1.75mm Midnight Blue, Reinforced 3D Printer Filament 1kg, Lightweight & High-Strength PETG-CF Filament, Strong & Easy to Print : Industrial & Scientific." Accessed: Feb. 19, 2025. [Online]. Available: <https://www.amazon.com/FLASHFORGE-Filament-Reinforced-Lightweight-High-Strength/dp/B0D7VSMCJS?th=1>
- [14] "Bambu Filament Technical Data Sheet V2.0". UW Makerspace. Accessed February 19, 2025.
- [15] "'CFN Medical Fiberglass Casting Tape | Medline." Accessed: Feb. 17, 2025. [Online]. Available: <https://www.medline.com/product/CFN-Medical-Fiberglass-Casting-Tape/Z05-PF140565?question=fiberglass&>
- [16] G. D. Smith, R. G. Hart, and T.-M. Tsai, "Fiberglass cast application," *The American Journal of Emergency Medicine*, vol. 23, no. 3, pp. 347–350, May 2005, doi: 10.1016/j.ajem.2005.02.031.
- [17] A. T. Berman and B. G. Parks, "A Comparison of the Mechanical Properties of Fiberglass Cast Materials and Their Clinical Relevance," *Journal of Orthopaedic Trauma*, vol. 4, no. 1, p. 85, Mar. 1990.
- [18] "Vacuum forming plastic - The ultimate guide," EuroPlas. Accessed: Feb. 19, 2025. [Online]. Available: <https://europlas.com.vn/en-US/blog-1/vacuum-forming-plastic-the-ultimate-guide>
- [19] International Committee of the Red Cross's Physical Rehabilitation Programme, *Manufacturing Guidelines: Ankle Foot Orthosis*. 2010.
- [20] D. Totah, I. Kovalenko, M. Saez, and K. Barton, "Manufacturing Choices for Ankle-Foot Orthoses: A Multi-objective Optimization," *Procedia CIRP*, vol. 65, pp. 145–150, Jan. 2017, doi: 10.1016/j.procir.2017.04.014.
- [21] "Polyethylene vs Polypropylene (PE vs PP): A Comprehensive Comparison Guide." Accessed: Feb. 17, 2025. [Online]. Available: <https://omnexus.specialchem.com/tech-library/article/polyethylene-versus-polypropylene>

Questions?

