

Kidney Clamp

for Laparoscopic Partial Nephrectomy

Jeff Hlinka, Shannon Hynes, Kayla Stankevitz Client: Dr. E Jason Abel, Department of Urology, UW- Madison Advisor: Mitch Tyler, Department of Biomedical Engineering, UW- Madison



Analysis and Results

· Solidworks analysis show displacement and stress distribution

THE T IL I I NVN // ///// N/

Abstract

Laparoscopic Partial Nephrectomy (LPN) surgery is indicated for removal of kidney tumors.. The current technique involves clamping the renal artery and vein. This procedure prevents all blood from entering the kidney tissue, causing the tissue to become ischemic which is risky for the patient. A new method of performing LPN surgery has been requested that would use a clamp to occlude blood flow in the kidney only to the portion of the organ being removed in an effort to maintain blood flow to the health tissue. An adjustable surgical clamp was designed and a prototype was manufactured using ABS. Preliminary testing was performed via solidworks analysis and it was determined that the device will not fail due to normal and maximum expected loading, 5 N and 20 N respectively. Further testing should be completed on a stainless steel prototype before conclusions can be drawn.

Motivation

- · More than 46,000 patients undergo nephrectomy surgery each year.
- · 32,000 new cases of renal cancer are diagnosed each year in the United States
- 1/3 of patients who undergo kidney resection will have a recurrence.
- · Clamping times of as little as 30 minutes have been shown to cause 10% loss in kidney function post-surgery2.
- · Market Need: a method to increase the effectiveness & safety of surgical treatment.

Background

The Kidney:

- · Filters around 20% of body's blood per minute.
- · Blood flow rate essential to maintain homeostatic functions (i.e. remove wastes).

Laparoscopic Partial Nephrectomy:

- · Removal of only diseased tissue from kidney.
- · Clamping needed to control blood loss and keep operating view clear

Current Method

Proposed Method Selectively occludes blood flow to the

Occludes blood supply to entire kidney portion of the kidney being removed causing cell death (global ischemia) by (regional ischemia) by clamping around the clamping at the source of blood flow. Can functional part of the kidney, or cause loss of kidney function if ischemic narenchyma for greater than 30 minutes.

Bulldog clama

Figure 1a (left) shows the current method of tumor removal with clamps applied at the blood source. Figure 1b(right) shows the proposed method, which employs a clamp around the kidney tissue1

Client

Our client Doctor E. Jason Abel at the University of Wisconsin - Madison hospital specializes in localized advanced kidney cancer. His philosophy to "provide maximal quality of life to patients by using minimally invasive approaches to cancer therapy," has prompted the idea for a new, laparoscopic tool to aid in partial kidney remova

Design Specifications

- · Must be able to fit through a 12 mm by 15 mm trocar · Provide enough strength to occlude renal blood flow
 - · Without causing harm to kidney
 - Adjustable force · Maintainable force
 - · Evenly distributed force
 - Maintain force for 5 30 minutes
 - Provide force 5 N or less³
- Reusable

Jaws

· . . .

•。][o

jaw in the 150° position

Yield Strength: 137.9 MPa

Young's modulus: 193.0 MPa

Ultimate Tensile Strength: 485 MPa

- · Able to operate with one hand
- Accommodate hand breadth ranging from 6.5 9.5 cm .
- . Not exceed cost of commercially available clamps (\$5000)
- Entire clamp should weigh less than 5 lbs.

· Made from spring stainless-steel (prototype from ABS)

Allows two conformations: 180° and 150°

Plunger deploys to lock clamp in 150° position

Figure 5 shows a Solidworks model of the jaws in the closed position. This position

allows the jaws to be flush and fit through the trocar. The enlargement shows the

plunger in the retracted position, to allow the top jaw piece to form a 180° angle

AISI 316 Annealed Stainless Steel

· 950x greater than our highest experimental stress

· Teeth provide grip to secure the organ

· Upper jaw has two pieces connected by hinge




Figure 2 shows the variation of force generated by current clamps on the market. The inlaid picture shows an example of the texture desired to stabilize the kidney

Final Design



Deployment Mechanism

 Jaws controlled via 2 wires · Wire to close top jaw relative to bottom jaw · Wire to retract plunger

Figure 3 shows the final design fabricated out of ABS using a 3D printer

· Shortened and released via handle

Arm 350 mm length

- 5 mm outer diameter
- · Familiar to surgeons
- · Connected to clamp

mechanism

Student Shop Employees

References

[1] Landman, Jaime. "Laparoscopic Partial Nephrectomy." Kidney Cancer Institute. Web. < http:// www.kidneycancerinstitute.com/PDFs/Acsculap-Nader-partialnephrectomy%20manuscript.pdf> [2] Peter L. Steinberg, Ravi Munver, and Reza Ghavamian. Journal of Endourology. September 2011 25(9): 1525-1530

[3] Farshad, Mehdi. "Material characterization of the pig kidney in relation with the biomechanical analysis of renal trauman." Journal of Biomechanics. 1999, 32: 417-425.

[4] Vahlensieck, W., Petrik, R. "Experiences with a new 'kidney holding and compression forceos," International Urology and Nephrology. 2(3). 253-260.

[5] Lee, H. Box, G., Abraham, J., Elchico, E., Panah, R. "Laboratory Evaluation of Laparoscopic Vascular Clamps Using a Load-Cell Device – Are All Clamps the Same?" Journal of Urology. 2008, 1267-1272

Figure 7: top images show force displacement of two top pieces modeled as one piece. Bottom image shows stress distribution from applied forces

Solidworks Analysis

Yield= 137.9 MPa	Normal Loading (5N)	Max Loading (20 N)	Extreme Loading (40 N)
Base Piece	30.66 MPa	122.63 MPa	243.76 MPa
Hinge Piece	16.34 MPa	65.35 MPa	130.70
End Piece	10.66 MPa	42.65 MPa	85.31 MPa
Plunger	13.42 MPa	53.67 MPa	107.34 MPa
Single Unit	15.33 MPa	61.30 MPa	122.61 MPa

Future Work

Stainless Steel Prototype:

- · Proof of concept complete
- · Fabricate stainless steel prototype Additional Function:
- · Need to design and fabricate a
- flexible arm
- ·Control arm flexibility
- · Incorporate RF ablation guide
- · Check for leak point pressure Test final prototype in pig lab.

· Verify occlusion force

occlude kidney

· Compare to pressure needed to

· Kidney ruptures at 20 N

· Ranges up to 5 N

Testing:

Acknowledgements

Dr. E Jason Abel, Department of Urology, UW- Madison Mitch Tyler, Department of Biomedical Engineering, UW- Madison

Heidi Ploeg, Department of Mechanical and Biomedical Engineering, UW- Madison John Puccinelli, Department of Biomedical Engineering, UW- Madison Sam Paulsen, BME Student, Roomate

• 90. 900 lbc •.031 cm Corrosion resistant •\$7.99/30 ft

allows the jaws to be parallel while clamping the kidney with 2 cm between the jaws. The enlargement shows the plunger in the deployed position to hold the top





Handle

· Ratchet to lock into place

through wire and bike brake

