

BME 301: Preliminary Product Design Specification

Team:	Mitchell	Josvai	Leader
	Kate	Eichstaedt	Communicator
	Avani	Lall	BWIG & BPAG
	Ashlee	Hart	BSAC

Advisor: Dr. Edward Bersu

Client: Dr. Amanda Doubleday

Date: February 12th, 2021

Function:

The function of the device is to aid in the retraction of the thyroid and parathyroid glands during thyroidectomies and related procedures. During surgery, endocrine surgeons must retract the thyroid gland medially in order to gain access to the recurrent laryngeal nerve and parathyroid glands. Depending on the procedure, they must then dissect the thyroid gland from vascular attachments, and possibly the parathyroid glands. Surgeons use stainless steel forceps with a piece of gauze clamped at the tip, referred to as a "peanut" to retract and hold the thyroid gland in place, without rupturing it. However, due to the single point of contact, the thyroid gland can often be too large to be held comfortably by this method. The device should be able to assist surgeons in retracting and holding the thyroid in place from multiple contact points.

Client Requirements

The client requires a surgical instrument to aid in the medially retraction of the thyroid and parathyroid glands during surgery. The device should have a single handle similar to a standard forceps, but with two prongs to retract the thyroid gland from multiple contact points. Each prong should be capable of clamping and holding a surgical peanut sponge, a small sponge used to reduce the forces on the thyroid. The handle of the device should have some sort of ratcheting system to be able to adjust the distance between the two prongs. Additionally, the ratchet should allow for the device to be held in place for a period of time without having to manually hold the clamps shut.

Design Requirements

1. Physical and Operational Characteristics

a. Performance Requirements:

The device should be able to assist in completing the tedious dissection of the thyroid without being traumatic. It should be reusable, auto-clampable and with blunt ends that act as clamps. It must function as one instrument that has two peanut forceps at the end opposite of the handle. The device must be capable of adjusting the width between the peanuts, and locking into the desired conformation. Finally, the device must be capable of withstanding all forces that are applied to it, both by the surgeon and the areas of the body it is acting on. The 95th percentile for human grip strength in right-handed men is around 500 N [1], and the device handle must be able to withstand this force. As the average adult's thyroid weighs between 20 and 30 grams [2], the forces applied by the thyroid are negligible in comparison to those applied by the surgeon.

b. *Safety:*

The only people allowed to operate with and use the device will be trained professionals, as to greatly reduce the risk of injury. In order to ensure the safety of the surgeons using the device, testing of the device's ability to endure forces of the hand, or ability to hold against the typical weight of a thyroid will be conducted. The device is not required to be permanently biocompatible, because it will only be in contact with the patient temporarily. However, the device must not be toxic, or susceptible to leaching of potentially harmful chemicals into the body. Finally, blunt edges and ends should be preferred over sharp edges, so as to avoid any unintended perforation or trauma caused by the device.

c. *Accuracy and Reliability:*

The device must reliably be capable of performing the task it is designed for. It must not puncture or cause trauma to the thyroid or other areas of the body when in contact. The ratcheting or latching mechanism must not jam or lock or when unintended, but also must lock when intended and avoid releasing until desired. The clamping mechanisms must be capable of holding onto a peanut for the length of surgery without risk of the peanut detaching and entering the body.

d. *Life in Service:*

The device should have a rather long life in service. The device should be able to be sterilized and reused by autoclaving. The "peanuts" or gauze pads that will be used on the ends of the device however are single use, but are not being designed by the team. If there are signs of device damage or material corrosion, the device should be replaced.

e. *Shelf Life:*

The device will be made out of surgical grade stainless steel. Due to the mechanical properties of stainless steel, the shelf life for the device will be rather long. The average lifetime of stainless steel products ranges from 15-25 years[3].

f. *Operating Environment:*

The device will be used in a surgical setting. The likely temperature that the device will be in is somewhere between room and body temperature, depending on the point in the procedure. This gives an operating temperature between 22° C and 37° C [4]. For pressure, the likely pressure the device will be experiencing is around 1 atmosphere [4]. The device will be used by surgical staff so it is important for the staff to receive adequate device training.

g. *Ergonomics:*

The device should be relatively simple to use by a trained user in an operational setting. The device should have a handle that is easy and comfortable to grip and is able to be held by one hand. The device should feature a ratchet that can vary the distance between the prongs so it does not have to be manually held to a certain distance. The device should not hinder the surgical staff during the operation, and ideally increases the ease of the procedure.

h. *Size:*

The current device being used is the Peanut Sponge Forceps and is approximately 8" in length [5]. Typical retractors used in surgery range from 8 to 12 cm. The device should be similar in length to these devices currently in use, so as to be easily adopted by

surgeons utilizing other methods of thyroid retraction. Thyroids are anywhere from 4-6 cm and the device should have the two prongs 2-3 cm apart on average, with an adjustable range spanning from 1-4 cm, so that it may be used on a variety of patients and thyroid sizes. Measurements will be taken on current surgical forceps and retractors to determine accurate dimensions.

i. *Weight:*

The weight of the device should be close to that of the weight of the forceps used with the peanut currently, or typical surgical forceps and retractors, at around .09 lbs [6]. A small increase of weight will be allowed due to the addition of the second prong, although the device should not be sufficiently heavy as to be difficult to operate by a surgeon. Measurements will be taken on current surgical forceps and retractors to determine a more accurate target weight.

j. *Materials:*

The device will be made out of stainless steel such as the current device and most modern surgical instruments are. Stainless steel provides greater durability because it is anti-bacterial, non-corrosive and rust-resistant. It is also autoclavable, which allows it to be sterilized quickly and repeatedly. The durable stainless steel construction means the device will last and remain dependable for medical use [5].

k. *Aesthetics, Appearance, and Finish:*

The medical device should have the appearance that of a typical surgical instrument. The device should have a highly polished, or mirror finish in order to prevent potential staining [7]. Other than this requirement, aesthetics are less critical to the design than other relevant criteria.

2. Production Characteristics

a. *Quantity:*

Only one device will be produced for the full project, but it will be reusable since it is made of stainless steel. In the future, if this kind of device is proved to be beneficial to the procedure more can be produced.

b. *Target Product Cost:*

The target cost of this device should be comparable to typical surgical forceps, although this cost varies greatly. Depending on the supplier and website, retail prices for surgical forceps and retractors can range from \$5.00 to around \$50 [8] from medical supply companies. For this design, the target cost of production will be between \$5.00 and \$10.00 per single thyroid retractor. Flexibility of cost will be taken into account to accommodate for the extra forceps incorporated into the device. Additionally, this is the final target cost of production, without development or prototyping taken into account.

3. Miscellaneous

a. *Standards and Specifications:*

If used during an actual operation, this product will need to be FDA approved because it falls under the category of a medical device which is protected under FDA regulation 21 CFR Part 807. If used during surgery, the device and the surgery must comply with CDC regulation regarding sterile procedures. During testing and clinical trials, the device must be tested under IRB regulations at the university level, and FDA regulation at the federal level.

b. *Customer:*

Our client has requested a two pronged, adjustable Peanut thyroid retractor. This device should ease common complications endured with just one prong, such as those associated with larger thyroids. Ultimately this device may suit other customers beyond our client, including other surgeons at UW and/or beyond.

c. *Patient-related concerns:*

To ease any patients' concerns, the device will be treated and used like any surgical device in the operating room. The device will be used by a trained professional, and cleaned thoroughly between uses.

d. *Competition:*

Currently, the Peanut is used but the one instrument is not enough traction and causes the thyroid to fold. To solve this problem, two Peanut Forceps are being used which is proving to also be problematic because it is difficult to maneuver both with one hand.

References

- [1] Stokes, J., 1976. NASA - MSFC-STD-512 - MAN/SYSTEM REQUIREMENTS FOR WEIGHTLESS ENVIRONMENTS | Engineering360. [online]
<https://msis.jsc.nasa.gov/sections/section04.htm>. Available at:
<<https://standards.globalspec.com/std/669461/MSFC-STD-512>> [Accessed 9 February 2021].
- [2] Pankow, B., Michalak, J. and McGee, M., 1985. *Adult Human Thyroid Weight*. [online] Pubmed.
Available at: <<https://pubmed.ncbi.nlm.nih.gov/4077513/>> [Accessed 8 February 2021].
- [3] “Stainless Steel for a sustainable future,” *Team Stainless*. [Online]. Available:
https://www.worldstainless.org/Files/ISSF/non-image-files/PDF/Team_Stainless_Stainless_Steel_for_a_Sustainable_Future.pdf. [Accessed: 10-Feb-2021].
- [4] Engineeringtoolbox.com. 2020. *STP - Standard Temperature And Pressure & NTP - Normal Temperature And Pressure*. [online] Available at:
<https://www.engineeringtoolbox.com/stp-standard-ntp-normal-air-d_772.html> [Accessed 14 September 2020].
- [5] “Peanut Sponge Forceps: Sklar Instruments 22-9480,” *quickmedical*. [Online]. Available:
<https://www.quickmedical.com/sklar-instruments-peanut-sponge-forceps.html>. [Accessed: 11-Feb-2021].
- [6] “ADC® Kelly Hemostatic Forceps, Straight, 5-1/2”L, Stainless Steel,” *Global Industrial*. [Online]. Available:
https://www.globalindustrial.com/p/medical-lab/medical-equipment/exam-room-supplies/kelly-hemostatic-forceps-straight-5-1-2-l-stainless-steel?infoParam.campaignId=T9F&gclid=Cj0KCQiApY6BBhCsARIsAOI_GjaErxyu_CezZTVpO3iKXoGy5DLCt760CsGWYqbcB1HmbmZV1jtzcEaApEOEALw_wcB. [Accessed: 10-Feb-2021].
- [7] “Metal Finishing for Surgical Instruments,” *RP Abrasives*, 15-Sep-2020. [Online]. Available:
<https://rpabrasives.com/industries/medical-surgical/>. [Accessed: 10-Feb-2021].
- [8] “MABIS Kelly Forceps, Medical Forceps, Locking Forceps, Silver, 5.5,” *Amazon*. [Online]. Available:
https://www.amazon.com/MABIS-Forceps-Medical-Locking-Silver/dp/B00EKQ7FY4/ref=sr_1_4?dchild=1&keywords=surgical+forceps&qid=1612984299&sr=8-4. [Accessed: 10-Feb-2021].